



1. BARRAQUER WIRE SPECULAM ADULT




ETS - 2301

2. PHACO WIRE SPECULAM ADULT & MEDIUM




ETS - 2302

3. WIRE SPECULAM SOLID BLADE




ETS - 2303

4. SELF RETAINING WIRE SPECULAM FOR SICS, CATRACT, PHACO & LASIK SURGERY




ETS - 2304

5. SELF RETAINING SOLID WIRE SPECULAM FOR WITHOUT ANEASTHESIA SURGERY



ETS - 2305




ETS - 2305A

6. CLARK EYE SPECULAM




ETS - 2306

9. PAEDIATRIC WIRE SPECULAM




ETS - 2309

7. MURDOCK EYE SPECULAM




ETS - 2307

10. ALPHONSA PAEDIATRIC WIRE SPECULAM




ETS - 2010

8. MULLER LACRAMAL RETRACTOR



ETS - 2308

11. IRRIGATING LASIK SPECULAM



ETS - 2311

12. DESMARES LID RETRACTOR
SIZE : 0,1 & 2




ETS - 2312

13. LID GUARD OR LID SPATULA




ETS - 2313

13. SCHIPENS ORBITAL RETRACTOR




ETS - 2314

15. MULLER'S RETRACTOR



ETS - 2315

16. STRABISMUS MUSCLE HOOK



ETS - 2316

17. JAMESON MUSCLE HOOK




ETS - 2317

18. B.P. HANDLE FLAT




ETS - 2318

19. B.P. HANDLE ROUND




ETS - 2319

20. TOOK'S CORNEAL KNIFE



ETS - 2320

21. FOREIGN BODY SPUD




ETS - 2321

21. Hockey Knife




ETS - 2321A

22. LENS EXPRESSOR VECTIS ROUND




ETS - 2322

23. LENS EXPRESSOR VECTICS SERRATED




ETS - 2323

25. ENUCLATION SPOON




ETS - 2324

24. EVESCREATION SCOOP



ETS - 2325

26. EVISCERATION SPOON




ETS - 2326

27. CHALAZION SCOOP




ETS - 2327

28. THERMAL HEAT CAUTERY BALL TYPE




ETS - 2328

29. PUNTUM DILATOR SHARP OR BLUNT TIP




ETS - 2329

**30. BOWMAN LACRAMAL PROBE SIZE:
0000, 000, 00, 0, 1, 2, 3 & 4**




ETS - 2330

31. PIG TAIL PROBE



ETS - 2331

32. HAMMER



ETS - 2332

33. PERIOSTAL ELEVATOR



ETS - 2333

34. LACRAMAL CHISEL



ETS - 2334

35. IRIS REPOSITOR DOUBLE END



ETS - 2335

36. HOCKEY KNIFE FOR CORNEAL DISSECTOR



ETS - 2336

37. LACRAMAL SAC DISSECTOR



ETS - 2337

38. FIXATION RING FOR PHACO SURGERY




ETS - 2338

**39. CASTROVELJO CORNEAL TREPHINE
SIZE : 4 MM TO 10 MM**



ETS - 2339

**40. DISPOSABLE CORNEAL TREPHINE
SIZE : 4 MM TO 10 MM**



ETS - 2340

41. CAPSULAR REXIS FORCEP UTRATTA




ETS - 2341

**42. MICRO CAPSULAR REXIS FORCEP UTRATTA
20G OR 23G**



ETS - 2342




**42. HALDIPURKHAR MICRO CAPSULAR
REXIS FORCEP**





ETS - 2342 A

**43. BIMUANUAL CANNULA FOR PHACO
HANPIECE 22 G IRRIGATION & 23G
ASPIRATION OR 21 G IRRIGATION & 22G
ASPIRATION**




RADIATION

ASPIRATION



ETS - 2343

**44. SILICON TUBE FOR SIMCOE
CANNULA & SLEEVE FOR
INSTRUMENTS**




ETS - 2344

45. I.O.L. DIALER




ETS - 2345

46. KUGLEN'S HOOK




ETS - 2346

47. IRIS REPOSITOR OR LENS SPATULA




ETS - 2347

48. NUCLUES 'Y' ROTATOR




ETS - 2348


49. CAPSULE POLISHER NIGHTINGALE




50. PHACO CHOPPER




51. NAGHARA SHARP CHOPPER




52. PHACO CHOPPER SHARP CURVED TIP



53. GREEN MUSCLE HOOK



54. KIRBY LENS EXPRESSOR



55. RETINAL DETACHMENT HOOK



ETS - 2355


56. DRESSING DRUM



Dressing Drum	6 x 6
Dressing Drum	9 x 9
Dressing Drum	11 x 9
Dressing Drum	14 x 12
Dressing Drum	15 x 12

ETS - 2356

57. TOWEL CLIP CROSS ACTION



ETS - 2357

58. BULL DOG CLAMP STRAIGHT




ETS - 2358

59. BULL DOG CLAMP CURVED




ETS - 2359

60. MOSQUITO ARTERY FORCEP STRAIGHT




ETS - 2360

61. MOSQUITO ARTERY FORCEP CURVED




ETS - 2361

62. SCISSOR STRAIGHT FOR DRAPE CUTTING




ETS - 2362

63. EYE SCISSOR CURVED



ETS - 2363

64. EYE LASH CUTTING SCISSOR CURVED



ETS - 2364

**65. CHALAZION FORCEP
SMALL, MEDIUM &
LARGE**



ETS - 2365

**66. CHALAZION FORCEP AYYER'S
PATTERN SMALL, MEDIUM &
LARGE**




ETS - 2366

**67. ENTROPIAM CLAMP
LEFT & RIGHT**




ETS - 2367

68. CALIPER STRAIGHT



ETS - 2368

**69. BONE PUNCH SIZE : 1MM,
1.5 MM, 2MM, 2.5 MM & 3MM**




ETS - 2369

70. BONE NIBBLING FORCEP




ETS - 2370

71. NASAL SPECULAM KILLIAN




ETS - 2371

72. BLADE BLEAKER PENCIL TYPE



ETS - 2372

**73. BLADE BREAKER
BARRAQUER TYPE**



ETS - 2373

**74. PIERSE HOSKIN FORCEP
STRAIGHT**




ETS - 2374

**75. PIERSE HOSKIN FORCEP
CURVED**




ETS - 2375

**76. BISHOP CORNEAL
FORCEP**




ETS - 2376

77. JAFFE CORNEAL TOOTHED FORCEP




ETS - 2377

78. LIMS CORNEAL FORCEP




ETS - 2378

79. CALLIBRI FORCEP




ETS - 2379


80. MCPHERSON CORNEAL FORCEP




81. ST. MARTIN CORNEAL FORCEP




82. SUPERIOR RECTUS FORCEP DASTOOR




83. EPILLATION FORCEP




84. IRIS FORCEP STANDARD MODEL




85. JEWEL TYPE TYING FORCEP




86. JAFFE TYING FORCEP STRAIGHT




87. JAFFE TYING FORCEP CURVED




88. MCPHERSON TYING FORCEP STRAIGHT



89. MCPHERSON TYING FORCEP ANGLED



90. MCPHERSON FORCEP



91. EYE DRESSING FORCEP



ETS - 2391

95. FIXATION FORCEP 2 X 3 TEETH




ETS - 2395

92. UTILITY FORCEP WITH CROSS SERRATED TIPS



ETS - 2392

96. KALT SUPERIOR RECTUS NEEDLE HOLDER



ETS - 2396

93. UTILITY FORCEP WITH VERTICAL GROOVED TIPS




ETS - 2393

94. FIXATION FORCEP 1 X 2 TEETH




ETS - 2394

97. BARRAQUER NEEDLE HOLDER FLAT




ETS - 2397

**98. BARRAQUER NEEDLE HOLDER
FINE FOR 10/0**




ETS - 2398

99. BARRAQUER NEEDLE HOLDER STRAIGHT




ETS - 2399

100. KELLY'S PUNCH




ETS - 2400

101. LENS HOLDING FORCEP KRATZ




ETS - 2401

102. LENS HOLDING FORCEP SHEPHARD OR IRIS CLAW CORCEP




ETS - 2402

103. LENS HOLDING FORCEP CLAYMAN




ETS - 2403

104. LENS HOLDING FORCEP BLADES




ETS - 2404

**105. MOODY'S FIXATION FORCEP FOR
SQUINT SURGERY**



ETS - 2405

110. TEFLON BLOCK



ETS - 2410

106. AKOSHI NUCLEUS PRE CHOPPER FOR PAHCO SURGERY



ETS - 2406

107. VITRECTOMY FORCEP WITH SMOOTH PLAIN TIP



ETS - 2407

108. VITRECTOMY FORCEP WITH SMOOTH SERRATED TIP




ETS - 2408

109. VITRECTOMY FORCEP WITH 1 X 2 TEETH




ETS - 2409

116. SCLERAL PLUG HOLDING FORCEP




ETS - 2416

117. CHARLES FLUTE NEEDLE CANNULA 20G




ETS - 2417

118. INFUSION CANNULA 20G - 2.5MM, 4MM, 6MM




ETS - 2418

**119. IRRIGATING VITRECTOMY LENS
PLANO - CONCAVE (45 DIAPTOR)**



ETS - 2419

120. IRRIGATING VITRECTOMY LENS BI-CONCAVE (90 DIAPTOR)




ETS - 2420

121. CORNEAL SCISSOR CASTROVEIJO SMALL BLADE




ETS - 2421

122. CORNEAL SCISSOR CASTROVEIJO MEDIUM BLADE




ETS - 2422

123. CORNEAL SCISSOR LEFT AND RIGHT




ETS - 2423

124. CAPSULE TAG CUTTING SCISSOR



ETS - 2424

125. BARRAQUER IRIS SCISSOR




ETS - 2425

126. VANNAS SCISSOR ANGLED FOR CAPSULATOMY



ETS - 2426

127. VANNAS SCISSOR STRAIGHT




ETS - 2427

128. VANNAS SCISSOR CURVED




ETS - 2428

129. WESTCOT CONJUNCTIVAL SCISSOR




ETS - 2429

130. SURGICAL BLADE NO. 11 & 15




ETS - 2430

**131. SIMCOE IRRIGATION & ASPIRATION
CANNULA SIZE : 18G, 20G, 21G,
22G & 23G LEFT HAND ASPIRATION**



ETS - 2431


**132. IRRIGATING VICTICS CANNULA - OVAL
SHAPE, SQUARE SHAPE & ROUND
SHAPE**



ETS - 2432




ETS - 2432 (A)



ETS - 2432 (B)


**133. SIMCOE IRRIGATION & ASPIRATION
CANNULA SIZE : 18G, 20G, 21G,
22G & 23G RIGHT HAND ASPIRATION**




ETS - 2433

134. CORTEX ASPIRATION CANNULA

ETS - 2434




**135. A.C. MAINTAINER CANNULA FOR
BLUE MENTHOL**




ETS - 2435

**136. 12'o CLOCK SIMCOE IRRIGATION &
ASPIRATION CANNULA LEFT & RIGHT**




ETS - 2436

**137. HYDRODISSECTION CANNULA
SIZE : 24G, 25G & 26G**



ETS - 2437

**138. AIR CANNULA 23G, 24G
25G, 26G & 27G**




ETS - 2438

139. 'J' HYDRODISSECTION CANNULA




ETS - 2439

141. DUCT LACRAMAL CANNULA



ETS - 2441

143. S.S. PERFORATTED TRAY WITH SINGLE SILICON MAT




ETS - 2443

145. INJECTION TRAY




ETS - 2445

140. HYDRODILENATION CANNULA




ETS - 2440

142. BANAJI LASIK CANNULA



ETS - 2442

144. S.S.PERFORATTED TRAY WITH DOUBLE SILICON MAT



ETS - 2444

145A. KIDNEY TRAY 6", 8" & 10"




Cheatle ETS - 2446




ETS - 2446 A

146. SPUN EYE DRAPE WITHOUT DRAINAGE
BAG SIZES : 60X45, 60X60 & 100 X 120



ETS - 2447

147 SPUN EYE DRAPE WITH DRAINAGE
BAG SIZES : 60X45, 60X60 & 100 X 120



ETS - 2448

148. TROLLEY DRAPE
SIZE 100 X 75




ETS - 2449

149. PATIENT'S LEGINS OR SHOE COVER




ETS - 2449A

150. POLY EYE DRAPE WITHOUT DRAINAGE SIZES 60X45, 60X60 & 100 X 120



ETS - 2450

151. POLY EYE DRAPE WITH DRAINAGE
SIZES 60X45CM, 60X60CM & 100 X 120CM



ETS - 2451

152. SURGEON GOWN



ETS - 2452

153. EYE PAD



ETS - 2453

154. PATIENT CAP



ETS - 2454

155. EYE SHEILD



ETS - 2455

156. SURGEON CAP & MASK



ETS - 2456

PHACO KNIVES

157. ENTRY KERATOME SIZES :
2.2 MM, 2.8 MM, 3.0 MM & 3.2 MM



ETS - 2457

158. TUNNEL BLADE SIZES :
2.2 MM & 2.6 MM




ETS - 2458
ETS-2458H

**159. LANCE TIP BLADE OR SIDE
PORT BLADE OR 15'o BLADE**




ETS - 2459
ETS-2459H

**160. SLIT KNIFE (EXTENSION
KERATOME) SIZES :**
5.1 MM OR 5.2 MM



ETS - 2460
ETS - 2460H



162. S.S. BLADE HANDLE




ETS - 2462

All required sizes are available, special sizes can also be made to order.
The above knives, are round stock knives & are manufactured from round wire stock.
This product has sharp edge as well as superb durability.
The below comparison states how these knives are better than flat stock knives.

NOTE : H-POLYCARBONATE HANDLE BLADES

Comparison



Flat Stock

Flat Stock knives are made on steel strips and from end these knives are wide at the time of cutting these blades can stretch the incision.




Round Stock

Round Stock knives are made on wire and at the time of from end these knives are wide at the time of cutting these blades can stretch the incision.




ETS - 2461
ETS-2461H

163. CARBON RAZER BLADE (RAJA BLADE)




ETS - 2463

164. BIPOLAR WET FIELD CAUTERY




ETS - 2464

165. BIPOLAR CABLE



ETS - 2465

166. BIPOLAR FORCEP CURVED




ETS - 2466

167. BIPOLAR FORCEP STRAIGHT




ETS - 2467

168. DISPOSABLE HYDRODISSECTION CANNULA




ETS - 2468

169. DISPOSABLE SIMCOE CANNULA 21G, 221G X 23G



ETS - 2469

170. DISPOSABLE SALINE CANNULA



ETS - 2470

171. TRIAL SET METAL RIM WITH
ILLUMINATION SIZE : 20 SPHERICAL X 6 CYLINDER



ETS - 2471

172. RETINOSCOPE MIRROR SINGLE




ETS - 2472

173. ARTIFICIAL EYE



ETS - 2473

174. NEAR VISION DRUM



ETS - 2474

176. CROSS CYCLINDER 0.26 & 0.5




ETS - 2476

175. RETINOSCOPE MIRROR
DOUBLE




ETS - 2475

177. TRIAL FRAME




ETS - 2477

178. ISHIHARA'S COLOUR VISION BOOK



ETS - 2478

179. 4X BINOMAC LOUPE WITH PLASTIC BAND OR ELASTIC BAND




ETS - 2479

180. A.R. PAPER ROLL




ETS - 2480

181. AUTOCLAVE INDICATOR TAPE ROLL OR SIGNAL LOCK




ETS - 2481

182. SCHOTZ TONOMETER



ETS - 2482

183. VISION DRUM WITH ILLUMINATION
(AVAILABLE BOTH 6 MTRS. OR
3 MTRS DISTANCE)



ETS - 2483

184. R.A.F. RULER




ETS - 2484

185. EXAMINATION SPECTACAL
(AVAILABLE 3X, 4X, 5X & 6X)




ETS - 2485

186. PINKY BALL



ETS - 2486

187. DIPLOBIA RED & GREEN GOGGLES




ETS - 2487

188. STREAK RETINOSCOPE



ETS - 2488

189. MADDOX SWING




ETS - 2489

190. BERGER LOUPE 4X



ETS - 2490

191. PRISSM BAR




ETS - 2491

192. LOOSE PRISSM SET




ETS - 2492

193. SUCTION TIPS (SET 4)




ETS - 2493

194. CHEATLE FORCEP




ETS - 2494

195. SPONGE HOLDING FORCEP




ETS - 2495

196. TILLEY'S FORCEP



ETS - 2496

197. SUTURE NEEDLE (18G, 20G)



ETS - 2497

CARE OF MICRO SURGICAL INSTRUMENTS

1. MANUAL CLEANING

- i. Immediately rinse the instruments with distilled water after use to remove blood, saline, tissue and other deposits which will stain the instruments./
- ii Always use a soft tooth brush with mild soap for cleaning.
- iii Wash the instruments thoroughly with distilled water and dry the instruments with a hot air blower or with a soft lint free cloth. (While using a cloth for drying be very careful as the tips can got damaged by contact)

2. ULTRASONIC CLEANING

- i. Ultrasonic cleaning is the best method of cleaning micro surgical instruments. Regular cleaning of the instruments with ultrasonic cleaner will increase the life of instruments.
- ii Strictly follow the instructions of the manufacturer of your ultrasonic cleaner & use only the recommended solution for cleaning.
- iii Do not allow the instruments to touch one another in the ultrasonic cleaner.
- iv Wash the instruments thoroughly with distilled water and dry the instruments with a hot air blower or with a soft lint free cloth. (While using a cloth for drying be very careful as the tips can get damaged by contact)

3. STERILIZATION

- i. Micro surgical instruments can be sterilized by the following methods a Autoclaving b. Dry heat c. Gas d. Dry chemical e.Boiling.
- ii Whichever method of sterilization is used, strictly follow the instructions of the manufacturer of the sterilizer
- iii Clean your instruments thoroughly before sterilizing, as sterilizing is not cleaning.
- iv. Instruments not required for surgery should not be sterilized, as sterilizing reduces the life of instruments.
- v. Regularly inspect the internal parts of your sterilizer for rust, as this rust can damage your instruments while sterilizing.
- vi. Do not allow the instruments to touch one another in the sterilizer.
- vii Do not allow the instruments to over heat, or be kept in the sterilizer for a long period of time, as this will damage the instruments.

4. INSPECTION

Inspect your microsurgical instruments regularly under a microscope to check the working parts & tips for damages. Do not use damaged instruments for surgery.


5. LUBRICATION

- i. Regularly lubricate the joints & locks of the instruments with a high grade lubricant.


6. STORING

- i. Clean the instruments thoroughly before storing
- ii Store the instruments in a dry place
- iii Always protect the tips of the instruments with soft silicone tubing before storage.
- iv Do not allow the instruments to move or touch one another when stored.
- v. Use trays with soft silicone mats for storing the instruments.
- vi All delicate and costly instruments, should be stored separately.


198. Hydrophilic Acrylic & Hydrophobic with Injector
360° Square Edge Aspheric Foldable IOL
2.2 mm x 2.4 mm (Preloaded also available)



Optic Size : 6.00 mm




Optic Size : 6.00 mm
 ETS - 2498




Optic Size : 6.00 mm

199. PMMA Capsular Tension Ring and Advanced Capsular Ring




Insertion Size : 9.00 mm
 ETS - 2499

200. Gloject Injector and Cartridge System




ETS - 2500


201. Intra Ocular Lenses PMMA Single Piece Posterior Chamber



Optic Size : 6.00 mm




Optic Size : 5.25 mm
 ETS - 2501




Optic Size : 5.50 mm

202. Intra Ocular Lenses PMMA Square Edge **203. Anterior Chamber IOL**



Optic Size : 6.00 mm
 ETS - 2502



Optic Size : 6.00 mm
 ETS - 2503